

GREEN IS CLEAN, USE CUROS

CUROS[®]

STERILE

R

Port Protector

Active Ingredient: 70% Isopropyl Alcohol (IPA) • Non-DEHP • Latex Free • Single Use Only ② • Rx Only

INSTRUCTIONS FOR USE

A **NEW** CuroS should be placed on the luer-activated needleless valve following completion of a valve access.



CuroS **Strip** Port Protector CuroS **Individual** Port Protector



1 Peel the foil seal from the CuroS cap.

2 Gently push and twist the CuroS clockwise onto the end of the valve.

Refer to current device labeling for indications for use, warnings and product details.

CUROS ALWAYS KEEPS PORTS PROTECTED & CLEAN

QUESTIONS? Contact Your Green Team ...



IVERA[™]
Solutions for a New Era in Infection Control[™]

Ivera Medical Corporation
3525 Del Mar Heights Road, Suite 430 • San Diego, CA 92130
888-861-8228 • www.iveramed.com
Patents Pending • Made in USA • ML-0028-9/11